

The background is a teal color with a complex, abstract pattern of lines. A dense network of thin black lines forms a grid-like structure that is slightly distorted, creating a sense of depth and movement. Overlaid on this are several thicker, curved black lines that sweep across the frame. A prominent, bright yellow line runs horizontally across the middle of the image, intersecting the black lines.

**NETZSCH**

Proven Excellence.

# UNLOCKING PYROLYSIS KINETICS

Practical Analysis with Kinetics Neo Software

## Erich NETZSCH GmbH & Co. Holding KG




### Analyzing & Testing

Devices for thermal analysis, determining thermophysical properties and rheology as well as fire testing systems



### Grinding & Dispersing

Extensive range of machines for wet and dry grinding, mixing, dispersing, homogenizing and separating



### Pumps & Systems

Always the right positive displacement pump for your application

Products and services for applications in the low and high temperature range from  $-260^{\circ}\text{C}$  to  $2800^{\circ}\text{C}$



## Thermal Analysis

Determination of dimensional and mass changes, phase transitions and enthalpies as a function of temperature



## Thermophysical Properties

Determination of thermal diffusivity and thermal conductivity, specific heat capacity and thermal expansion coefficients



## Adiabatic Calorimetry

Analysis of decomposition processes and reaction processes with regard to temperatures, released heat quantities and pressure curves



## Rheology

Measurement of rheological properties of non-Newtonian fluids and soft solids - from formulation to product use



## Fire Testing Devices

Determination of the fire behavior of products in the automotive, construction, electronics and polymer sectors; classification into "European fire classes"



**>150M €**

Sales



**>700**

Employees worldwide



**>10**

Laboratories worldwide



**6**

International  
production sites



**> 50**

Sales and service  
locations worldwide

# NETZSCH Analyzing & Testing - worldwide presence

# NETZSCH

Proven Excellence.



- 1. Pyrolysis and current challenges**
- 2. Kinetics in Solids and Kinetics Neo software**
- 3. Review Example Study of the Olive Stone Biomass Pyrolysis**
- 4. Software Demonstration**



# Pyrolysis and current challenges

- Pyrolysis is the thermal **Decomposition reaction** of organic compounds in an inert atmosphere.
- During pyrolysis, solid, liquid or gaseous products can be generated. If gases are released from the sample during pyrolysis, the changes in mass can be detected by [TGA \(thermogravimetry\)](#).

## Materials for pyrolysis

- Biomass
- Plastics
- Rubber

- Thermogravimetric Analysis (TGA)
- Stages of Decomposition (Pyrolysis)

## Kinetics Analysis


## Products from Pyrolysis

- Gases
- Liquids
- Solid Residue (Charcoal)

- The poor choices of the kinetic methods and experimental conditions
- kinetic complexities of multi-step pyrolysis processes



# Kinetics in Solids



Kinetics is the study of the dependents of a chemical reaction rate on time and temperature.

General equation

$$\frac{d\alpha}{dt} = k(T) \cdot f(\alpha)$$

Arrhenius dependence

$$k(T) = A \cdot \exp\left(\frac{-Ea}{RT}\right)$$

$$\frac{d\alpha}{dt} = A \cdot \exp\left(\frac{-Ea}{RT}\right) \cdot f(\alpha)$$

$\alpha$  : conversion


A: pre-exponential factor [1/s]

Ea: activation energy [kJ/mol]

R: gas constant 8.31 [J/(mol.K)]


T: absolute temperature [K]

$$T[\text{K}] = T[^\circ\text{C}] + 273.15$$



- **Mass sample 10 mg**
- At least 3 measurements at **different heating rates with  $R^2$  value greater than 0.995**
- Heating rates of **1, 2, 5, 10, and 20 K/min** are the most commonly reported in the literature for pyrolysis.
- **Reproducibility**






## TG 309 Libra Classic



Thermogravimetry: mass change is measured during heating

# Conversion $\alpha(t)$ for TGA data

**TG 309 Libra Classic**



$$\alpha = \frac{m_0 - m_t}{m_0 - m_\infty}$$


$m_0$ : initial mass  
 $m_t$ : mass at time  $t$   
 $m_\infty$ : final mass

- $\alpha(t)$  : Conversion is the ratio of the partial mass loss at given time point to the total mass loss at the final time point



# Kinetics Neo Software

- **Fitting Experimental Data by using Mathematical models:** Mathematical models enable the fitting of experimental data to theoretical equations, which helps in determining kinetic parameters like **Activation energy**, **Pre-exponential factor**,  **$f(\alpha)$**  and **Coefficient of determination  $R^2$** .
- **Kinetics Neo software** can simulate and **predict** the behavior of chemical processes under various temperature conditions, such as isothermal, dynamic, and adiabatic profiles.
- **Kinetics Neo Software** allows users to **optimize** temperature programs.



# Supported Measurement Types



ARC / MMC



Rheology



DMA



DIL



- NETZSCH and non-NETZSCH data
- ASCII text files like TXT or CSV
- Data Columns ( Temperature, Time and Signal )



DSC



DEA



TGA

**Data import**

**Single-point model-free**

**Iso-conversional model-free**

**Model-based kinetics**

**Statistical comparison**

**Prediction**

**Rate control**

**Model design**

**Curve fit**

**Manual optimization**

**Kinetics parameters**

**Automatic optimization**

NETZSCH Kinetics Neo - Glucose\_Analysis.kinx

Model Based

Mass / %

Temperature / °C

Legend

- 10.3 K/min
- Fit
- 5.1 K/min
- Fit
- 2.0 K/min
- Fit
- 1.0 K/min
- Fit

Model t:

Describe: Fn, Fn

Optimize Fit To: Signal

Reaction Steps

A → B   Cn → d →

B → C   Fn → d →

B → D   Fn → d →

Equation

$$d[b] \rightarrow d[d] = \text{PreExp} \cdot b^n \cdot \text{Exp}[-\text{ActivationEnergy}/(RT)]$$

parameters	Show Range	Value
ActivationEnergy		170.965
Log(PreExp)		13.388
ReactOrder n		1.940
Contribution		0.738



Recalculate   Optimize

Mass   Show Range

GLUCOS10   Value

4.293

# Multi-step Model-fitting: Connection Between Steps



## Independent


Mixture of non-interacting components

## Consecutive

Reaction in the single component

## Competing


Product depends on heating rate





**Independent:** reactant and product are not involved in other steps

**Consecutive:** product of one step is a reactant of another step

**Competing:** involve the same reactant






**Coefficient of determination  $R^2$  or R-squared** : The value  $R^2$  expresses the deviations between the calculated model curves and the experimental data compared to the deviations of the experimental data from their average value.



Without  $R^2$ , kinetic analysis lacks quantitative validation and is more prone to error.  
 $R^2$  confirms model accuracy and ensures reliable kinetic predictions.

**Curve fitting visualization:** Visualization helps validate how well a model fits the data and reveals errors like outliers or poor model choices



R<sup>2</sup> is a helpful but it doesn't tell the whole story. Without curve fitting plots, you're "fitting blind"

# Step-by-Step Guide to Kinetic Analysis

Unlocking the Pyrolysis of Olive Stone Biomass: TGA Analysis and the Appropriate Kinetic Approach




# How to Determine Thermal Stability of an Olive Stone With TGA?

## TG 309 Libra Classic



3.3% mass loss moisture evaporation




TGA measurement on olive stone at different heating rates; solid lines: TGA, dashed lines: DTG


Instrument	NETZSCH TG 309 Classic
Crucible	Al <sub>2</sub> O <sub>3</sub> , open
Sample mass	9.65 mg to 9.85 mg
Temperature range	25°C to 1000°C
Atmosphere	Nitrogen (40 ml/min), switch to synthetic air (40 ml/min) at 900°C
Heating rates	2.5 K/min, 5 K/min, 7.5 K/min, 10 K/min, 15 K/min, 20 K/min,

TGA measurement

→ decomposition temperature






# Measurements at Different Heating Rates



Decomposition of olive stone to 900°C at different heating rates, measured TGA data

How to Determine the number of steps?







Decomposition of olive stone to 900°C at different heating rates, measured TGA data

Conversion rate of the measurement at 2.5 K/min to 700°C. One peak and 4 shoulders indicate a 5-step decomposition process


$$\frac{d\alpha}{dt} = k(T) \cdot f(\alpha)$$

Source Data



- Bouzbib, M., Moukhina, E., & Schmölder. *Unlocking the Pyrolysis of Olive Stone Biomass: TGA Analysis and the Appropriate Kinetic Approach.* [www.netzsch.com](http://www.netzsch.com)

# Kinetic evaluation of the decomposition of olive stone



Kinetic parameters of the thermal degradation of olive stone

Reaction step	A → B Fn <sup>1</sup>	B → C Fn <sup>1</sup>	C → D Fn <sup>1</sup>	D → E Fn <sup>1</sup>	F → G DFn <sup>2</sup>
Activation energy [kJ/mol]	151.824	165.479	194.592	206.720	179.468
Log (Pre-Exp) Log (1/s)	14.083	13.792	15.116	15.286	12.093
Reaction order	1.832	2.732	1.039	1.466	6.304
Contribution	0.061	0.336	0.313	0.073	0.217
Coefficient of determination	0.999				


<sup>1</sup>Fn: Reaction of n<sup>th</sup> order

<sup>2</sup>DFn: One-dimensional diffusion of n<sup>th</sup> order


100 200 300 400 500 600 700 800 900  
Temperature / °C

Kinetic evaluation of the decomposition of olive stone. Rhombus lines: measured curves; solid lines: calculated curves based on a five-step reaction

## Conversion 5 %/min




## Conversion 1 %/min




Mass-loss control for process optimization at the same temperature program: solid line (mass loss %) and dashed line (temperature)

Mass-loss control for process optimization at the same temperature program: solid line (mass loss %) and dashed line (temperature)



Prediction of mass change over time under different isothermal conditions



Prediction of conversion over time under different isothermal conditions

Properties

**MultipleStep Prediction**

Method / Model  
p; Fn Fn Fn Fn DFn


Start T., °C	End T., °C	Hea Rat K/min	Time min		
10.0	460.0	5.00	90.0	+	×
460.0	900.0	2.00	220.0	+	×
900.0	900.0	0.00	150.0	+	×

Add Step Import Export

Calculate

Show additional curves

Temperature program



Multiple step prediction of mass change over time under different temperature conditions




# Software Demonstration

# Step-by-Step Recap

**NETZSCH**

Proven Excellence.

**TG 309 Libra Classic**



At least **3** measurements at different heating rates.

Materials for pyrolysis

Measured data TGA

## Kinetics Analysis

- Import data
- Create the most suitable model
  - Reaction type for each step
    - nth order...

Validation of kinetic models

**Kinetics Triplet, Prediction and Process Optimization**



**ICTAC:** *International  
Confederation for  
Thermal Analysis and  
Calorimetry*



- Model free analysis
- Multi-step model-fitting (model based)
- Diffusion control for curing
- Crystallization kinetics
- Kamal model for curing
- Deconvolution analysis (sum of peaks)



Thermochimica Acta  
Volume 689, July 2020, 178597



Review

## ICTAC Kinetics Committee recommendations for analysis of multi-step kinetics

Sergey Vyazovkin <sup>a</sup>, Alan K. Burnham <sup>b</sup>, Loic Favregeon <sup>c</sup>, Nobuyoshi Koga <sup>d</sup>,  
Elena Moukhina <sup>e</sup>, Luis A. Pérez-Maqueda <sup>f</sup>, Nicolas Sbirrazzuoli <sup>g</sup>

<sup>a</sup> Department of Chemistry, University of Alabama at Birmingham, 901 S. 14th Street, Birmingham, AL, 35294, USA

<sup>b</sup> Alan Burnham Consultant, 4221 Findlay Way, Livermore, CA, 94550, USA

<sup>c</sup> Mines Saint-Etienne, University of Lyon, CNRS, UMR 5307 LGF, Centre SPIN, F-42023 Saint-Etienne, France

<sup>d</sup> Department of Science Education, Graduate School of Education, Hiroshima University, 1-1-1 Kagamiyama, Higashi-Hiroshima 739-8524, Japan

<sup>e</sup> NETZSCH-Gerätebau GmbH, Wittelsbacherstrasse 42, Selb 95100, Germany

<sup>f</sup> Instituto de Ciencia de Materiales de Sevilla, C.S.I.C.–Universidad de Sevilla, C. Américo Vespucio No. 49, 41092 Sevilla, Spain

<sup>g</sup> University Côte d'Azur, Institute of Chemistry of Nice, UMR CNRS 7272, 06100 Nice, France


Received 18 March 2020, Accepted 19 March 2020, Available online 16 May 2020, Version of Record 5 June 2020.

■ Go to: <https://kinetics.netzsch.com>

■ Trial Version 30 days



[kinetics.neo@netzsch.com](mailto:kinetics.neo@netzsch.com)



## Kinetics Neo

Software for Kinetic Analysis and Simulation of Thermoanalytical Data for Chemical Processes

Kinetics Neo software fully supports 'ICTAC Kinetics Committee recommendations for analysis of multi-step kinetics'.

### News

**October, 23, 2024 Webinar Ceramics Sintering: Kinetics, Simulation and Process Optimization Using the Kinetics Neo and Termica Neo Software.**

In this webinar, we will present the typical solution steps sintering optimization of different ceramic materials in order to get the highest quality at lowest costs. They include kinetic modelling of the process by the NETZSCH Kinetics Neo software and then the simulation of this process for the user's geometry by the [Termica Neo](#) software. This helps understand your process and saves a lot of time and efforts compared to the way of trial-and-error.  
[Register for Webinar](#)

**September 25-27, 2024.** Kinetics Neo and Termica Neo will be presented on 50<sup>th</sup> GEFTA annual conference 2024 in Gießen, Germany is 2 talks:

- Elena Moukhina, Jan Hanss. Kinetic Modeling of Metal Reduction at Different Temperature Conditions and Hydrogen Concentrations

You can rely on NETZSCH.

**NETZSCH**  
Proven Excellence.

# Thank you for your attention!

**Dr. Mohammed Bouzbib**

Chemist

For further questions please contact

[kinetics.neo@netsch.com](mailto:kinetics.neo@netsch.com)  
[webinar\\_ngb@netsch.com](mailto:webinar_ngb@netsch.com)

kinetics.netsch.com

**NETZSCH**

Proven Excellence.

150

EST  
1873